
Civière d'évacuation / Rescue stretcher / Evakuierungstrage
/ Camilla de evacuación / Barella di evacuazione / Evacuatie
brancard / Nosze ewakuacyjne / Maca de evacuação


FA 70 018 00

KRATOS SAFETY


689 Chemin du Buclay
38540 Heyrieux - FRANCE

Tel : +33 (0)4 72 48 78 27
Fax : +33 (0)4 72 48 58 32


www.kratossafety.com info@kratossafety.com

1


2


3

4


5



6

Déplacement vertical / Vertical movement / Vertikale Beförderung / Desplazamiento vertical / Spostamento in verticale / Verticale verplaatsing / Przemieszczanie w pionie / Transporte vertical:

Déplacement horizontal / Horizontal movement / Horizontale Beförderung / Desplazamiento horizontal / Spostamento in orizzontale / Horizontale verplaatsing / Przemieszczanie w poziomie / Transporte horizontal:

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This guide must be translated by the dealer into the language of the country of use (except if the translation is provided by the manufacturer).

For your safety, comply strictly with the instructions for use, verification, maintenance and storage.

KRATOS SAFETY cannot be held liable for any direct or indirect accident occurring following use other than that provided for in these instructions; do not use this equipment beyond its limits!

USE AND PRECAUTIONS:

The KRATOS SAFETY stretcher is a flexible stretcher designed for a variety of uses, however it is mostly intended for rescue operations in restrictive environments, confined spaces and for vertical rescues. The stretcher comes with a carry bag containing all the accessories required for rescue operations.

Technical specifications:

Maximum weight of the victim: 150 kg

Weight of the stretcher: 5.5 kg

Weight with bag and accessories: 8 kg

Stretcher material: Polyethylene

Strap material: Polyamide

Temperature for use: - 50 °C up to +70°C

The bag contains:

- ▶ The stretcher, with 4 black horizontal carrying handles (2 on each side), 4 black safety straps on the side with sewn-on loops and brass grommets, 2 black safety straps at the bottom and 1 red safety strap at the top
- ▶ 2 connectors for attaching the stretcher to the evacuation system, and for any additional stretcher attachments
- ▶ 2 yellow straps for moving the stretcher horizontally
- ▶ 4 black/white straps for creating 4 additional handles
- ▶ 2 black/white straps for additional safety
- ▶ 1 black strap with built-in handles for pulling the stretcher on the ground.
- ▶ 1 9-m long Kernmantle rope for moving the stretcher vertically.

Warning! NEVER suspend the stretcher using the grommets.

Use the cords and straps provided in the original bag and according to their intended purpose.

This product must only be used by trained and capable individuals.

Positioning the victim:

The stretcher is not intended to immobilise the spine. If a spine injury is suspected, protect the patient using a spinal immobilisation device before moving him/her onto the stretcher.

- 1/ Take the stretcher and accessories out of the bag.
- 2/ Unroll the stretcher, and if necessary roll in the opposite direction so that it stays flat once placed on the ground. Put the accessories to one side and position the stretcher as shown below:
- 3/ Release the 4 black safety straps on the side, the two black safety straps at the bottom and the red strap at the top to make room for the victim.
- 4/ Slide the victim onto the stretcher:
 - **Feet first:** cross the victim's arms, start by sliding the top part of the stretcher (the side with the red strap) under the victim's feet then slide the victim along the stretcher lengthwise, pulling by the legs if possible.
 - **Head first:** cross the victim's arms, start by sliding the bottom part of the stretcher (the side with the black safety straps) under the victim's head then slide the victim along the stretcher lengthwise, pulling by the shoulders if possible.
- 5/ Attach the straps working from the top towards the bottom, buckle the 4 side safety straps then the 2 black safety straps at the bottom. These must be looped through the grommets near the last black side safety strap. If necessary, adjust the tightness of the black safety straps; the victim must not have any freedom of movement, make sure there are no compression points on the victim.

Make sure that the victim's feet are fully on the stretcher.

Do not cross the side safety straps, they must be horizontal at all times.

Carrying the stretcher with the victim:

6/ Movement using human strength:

- A victim placed on a stretcher can be carried by 4 rescuers.

If necessary, additional handles can be added. To do this, use the black/white straps, fold them in half, make a knot at the end and pass them through the unused grommets on the side.

This will allow the stretcher to be carried by 8 rescuers if necessary.

- A victim placed on a stretcher can also be pulled over the ground by 3 rescuers:

1 in front using the red strap, and 2 on the sides using the black strap with built-in handles.

7/ Horizontal movement:

When moving the stretcher horizontally, the 2 yellow straps provided for this purpose must be used. These straps must pass under the stretcher and come out through the angled slots.

Warning! The yellow straps are not the same length, the shortest one must be placed at the front part of the stretcher and the longest at the back.

These straps must be made using the provided connector, and attached to the rescue system (or a helicopter-hoisting device). Check that the connector is tight before starting any rescue operation.

When the stretcher is moved in this way, the red strap at the top must be attached to the second side safety strap to create a "safety hood" for the victim's head.

8/ Vertical movement:

When moving the stretcher vertically, the 9-m Kernmantle rope provided for this purpose must be used. The rope must be passed through the unused side grommets.

When the stretcher is moved in this way, the red strap at the top must be attached to the second side safety strap to create a "safety hood" for the victim's head.

Re-packing the stretcher:

Check the integrity of all the stretcher's components. In particular, check that there are no marks, cracks, signs of wear, cuts, burns, frayed threads, etc. Check the condition of the cord ends and the connectors.

If you are in any doubt as to the condition of any of the components in the kit, do not use it and/or return it to the manufacturer or a competent person authorised by the manufacturer.

If all the components are in perfect condition, re-pack the stretcher the way it came.

These instructions for use must be kept with the product, inside the pouch provided and safe from any substance or liquid that could compromise perfect legibility.

MAINTENANCE AND STORAGE: (These instructions must be strictly observed)

During transportation, keep the equipment away from any cutting edges and in its packaging. Clean with soapy water, dry with a dry cloth and hang it in a naturally ventilated area away from any direct fire or heat source, as well as any elements which became wet during their usage. Bleaches and detergents are strictly prohibited. The equipment must be stored in a temperate, dry and ventilated area in its packaging **protected from sunlight** (which could affect the synthetic fibres), heat and chemical products.

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